

# TurnAid Pullsheet



Material: 100% Cotton

## Use:

TurnAid Pullsheet is used together with TurnAid Slide sheet for TurnAid turning system T2 and T3, and is mounted on the 2 rollers of the system.

## Placement of Turning Sheet:

Turn on the TurnAid.

Fold both sides of the system so both motor rolls will turn upwards.

Turn each motor roll so the Velcro tape appears on the top side of the roll.

Place the turning sheet in the middle of the bed, with the Velcro tapes on each side of the sheet along the rolls. REMEMBER to place the sheet with its bottom side downwards, so its Velcro tapes also turn downwards.



Starting with the upper sheet corner, attach the Velcro tape to the roll in one sliding movement while your hand straightens the sheet all the way down to the bottom end.



The sheet must turn across the top side of the roller, to ensure a correct direction of the sheet when you operate the rolls. A wrong placement may cause the patient's arms, fingers or others to be trapped between the sheet and the roll. The control only provides a low power to the motors in the tightening direction and they will stop in case there is any overload in this direction.



Tighten in approx.  $\frac{3}{4}$  of the sheet on the motor roll by pressing  or  on the hand control.

Attach the turning sheet to the Velcro tape on the opposite motor roll.



Tighten the turning sheet on the motor roll by pressing  or  until it is straight on the mattress with its dark stitching in the middle of the mattress.

RECOMMENDATION: If possible, keep the sheet's dark stitching in the middle of the mattress. This ensures the same amount of sheet on both sides.

## Cleaning:



80° Machine wash. Shrinks app. 4%



Tumble dry at normal temperature



Ironing at middle temperature